

Classical Isothermal Microcalorimetry (for radionuclidic standardizations)

R. Collé

research chemist (retired)

formerly of NIST (1974 – 2003)



VERMI

Virtual European Radionuclide Metrology Institute

Young Researchers Workshop

1–5 December 2003 – Paris

... I have not much doubt but that the two standards will be found in very good agreement, but it will be a devil of a mess if they are not.

... I think I can compare two nearly equal standards [to] an accuracy of 1 in 1000. I suppose, however, we shall not worry if the agreement of the two standards is within 1 in 300 or 400.

Ernest Rutherford (University of Manchester)

letter to Bertram Boltwood (Yale University)

18 march 1912

Some quick history

- 1896 radioactivity discovered
- 1903-4 first primary (“absolute”) measurements of radioactivity (?)
- P. Curie & Labord (1903) – first calorimetry of radioactivity
- P. Curie & Dewar (1904) – cryogenic calorimetry at liquid water & liquid hydrogen temperatures
- Rutherford & Barnes (1903-4) dual cell & differential calorimetry
- 1914 – 45+ calorimetry was a principal tool for international intercomparison of standards (radium)
- 1925 Ellis & Wooster – confirmed need for neutrino
- > 1950 calorimetry of “special nuclear materials” (power & bombs)
- Now ! milli-K calorimetry / quantum bolometry

But you can still do old-fashioned calorimetry

Some references – classical calorimetry


- O.E. Meyers, *Nucleonics* **5** (11), 37 (1949) -- exhaustive review & bibliography of calorimetric measurements of radioactivity up to 1949.
- W.B. Mann, *Encyclopedic Dictionary of Physics*, Pergamon, Oxford (1962) – updates Meyer's review to 1958.
- H.L. Callandar, *Proc. Phys. Soc.* **23**, 1 (1910) – Peltier effect “radiobalance” for comparison of international standards
- W.B. Mann, *Nucl. Instr. Meth.* **112**, 273 (1973); *Appl. Radiat. Isot.* **46**, 185 (1995) – the NIST radiobalance in use from 1953 to 1980s
- ➔ • H. Ramthun, *Nucl. Instr. Meth* **112**, 265 (1973) – GREAT survey of different types of “classical” calorimeters
- R. Collé, *Appl. Radiat. Isot.* **56**, 223 (2001). Classical isothermal microcalorimetry at cryogenic (8 K) temperatures
- R. Collé, coming soon ... use of a commercial microcalorimeter

Radionuclidic Microcalorimetry

WHY ?

Need to standardize GBq sources

Non-destructive (?)



dual-compensated cryogenic (8 K)

microcalorimeter

Appl. Radiat. Isot. 56, 223-230 (2002)



dual-cell near-isothermal (heat flow)
microcalorimeter (303 K)

Recent calorimetric-based standardizations 2000-2003

- ✚ verified extant calibration factors for (i) *Radiance* ^{32}P “hot-wall” angioplasty balloons and (ii) *Novoste* old-style, ceramic-cored, ^{90}Sr - ^{90}Y intravascular seeds
- ✚ performed primary standardization for *Novoste*, new-generation, aluminum-cored ^{90}Sr seeds to establish calibration factors
- ✚ primary standardization for ^{103}Pd (and for calibration of *Theragenics* prostate seeds)

Basic relationship between

Rate of energy (heat) input , or power P ,

and

Activity A


is
$$dH/dt = P = A \hat{E}$$

\hat{E} = average energy per decay

$^3\text{H} / ^{55}\text{Fe}$	0.9 $\mu\text{W}\cdot\text{GBq}^{-1}$
$^{103}\text{Pd} / ^{125}\text{I}$	9.
^{32}P	111.
$^{90}\text{Sr}-^{90}\text{Y}$	181.
^{226}Ra	4338.

Assumes absorb & measure
ALL ionizing radiation (no
losses)


And no “heat defect” effects
(i.e., no chemistry)



$$P = dH/dt = -k (T - T_e)$$




constant




$$T_S = T_R$$


$$P_S = dH/dt = P_R$$




CSC “Isothermal Microcalorimeter (IMC)”














**port assemblies -- source
(absorbers) holders & cells**

$$P = V^2/R$$






$$P = V^2 / R$$




$$\text{mean}(P) = 1/2 \text{ mean}(P_{+V}) + 1/2 \text{ mean}(P_{-V})$$

$$\begin{aligned} \text{var}(P) = & 1/2 \text{ var}(P_{+V}) + 1/2 \text{ var}(P_{-V}) + \text{covar}(P_{+V}, P_{-V}) \\ & + \text{autocorr}(P_{+V}) + \text{autocorr}(P_{-V}) \end{aligned}$$



typical calibration
data set




Filliben normality test
for calibration data
(Novoste seeds)

COMBUSTED BALLOON "A" ^{32}P

RUN #35

HAD
K=5882
±10
FROM
RUN
34A
p.63

P
K=5883
MW





IMC Lab Book II, p.79

combusted ^{32}P balloon "A"

RUN #36

$\frac{P_{5883}}{\text{MW}}$



IMC Lab Book II, p.97

uncombusted ^{32}P balloon "C"

Standardization scheme for Novoste ^{90}Sr new-style seed

19 seeds in 3 batches

$\left\{ \begin{array}{l} \text{X1 - X8} \\ \text{Y1 - Y8} \\ \text{Z1 - Z3} \end{array} \right\}$

ION CURRENT MEAS.

On all 19

PIC "A" $\Gamma_I = I_i / I_{\text{RRS}}$
 Vinten (NPL) I_i

CALORIMETRY $P = \bar{E} A$

ESTABLISH CALIBRATION FACTOR

$K = \Gamma_i / A_i$

16-Seed Train Returned (X + Y) s.

Meas. In batches

All 8 Xs [2] $P_X \rightarrow A_X$
 All 8 Ys [2] $P_Y \rightarrow A_Y$
 All 16 (X + Y)s [3] $P_{XY} \rightarrow A_{XY}$
 (Z1 + Z2) [2] $P_Z \rightarrow A_Z$

w/ Γ_I assign individual seed activity
 $A_{Xi} = (\Gamma_{Xi} / \sum_i \Gamma_{Xi}) A_X$

CHEMICAL DIGESTION / LS ASSAY

$A_{Z3} = \Gamma_{Z3} / K_i$

COMPARE

Z3

DIGEST
RINSE
DILUTE

M

DILUTE

S

DISPENSE


ampoules

LS vials


FROM MASSIC ACTIVITY ASSIGN A_{Z3}

Residual activity on tools by LS

2 scintillants 3 spectrometers
 ASSAY w/ 3H efficiency tracing

Novoste ^{90}Sr new seeds (Z1+Z2)


Radiance ^{32}P balloon "A" (combusted)



Get P_B and P_0 from "fit"

$$P = P_B + P_0 \exp(-\lambda t)$$

Novoste ^{90}Sr - ^{90}Y new seeds (16)



requires baseline P_B
measurement

Standardization results

source batch	number in batch	power (micowatt)	equivalent activity (MBq)
all (Xs + Ys)	16	342.70 342.30 344.41	1893.2 1891.0 1902.6
all Xs	8	172.45 172.76	952.7 954.4
all Ys	8	169.84 170.75	938.3 943.3
(Z1 + Z2)	2	57.81 58.27 58.09	319.4 321.9 320.9

From ion current measurements on individual seeds

X1 -X2 ranged from 114.4 MBq to 121.7 MBq
 Y1 -Y2 ranged from 113.9 MBq to 119.9 MBq

} $\pm 1.6 \% (k = 2)$

Z1 -Z3 ranged from 156.3 MBq to 164.3 MBq $\pm 3 \% (k = 2)$

Destructive assay (by LS) on Z3 agreed with calorimetry
initially by - 1.1 % (now with n = 5 determinations to 0.4 %)

Novoste new-style (Z batch) ^{90}Sr seeds

at power = 60 microwatt

number of determinations	agreement with LS-based destructive assay
1	1.8 % *
3	0.7 %
7	0.4 %

* typical 1-2 % accuracy for a single
determination at this power level

FINDINGS

**Got 1 to 2 % agreement w/ extant IC calibrations for both
Radiance ^{32}P balloons
Novoste ^{90}Sr - ^{90}Y seeds (old)**

**Calibration factors (15 determinations) for Novoste seeds (new)
has s.d.m $< 0.1\%$ and $1/2\text{range} = 0.5 \%$**


**Better than 0.5 % agreement between LS-based digestive assay
for new-style Novoste ^{90}Sr - ^{90}Y seeds**

**Replicate measurement uncertainty is about 0.5 % to 1 % or so
if can get P_B by fit with decay
or if one has sufficient replications to get ΔP
(with little decay)**

need primary standardization

+ calibration factors


+ transfer standards




Brachytherapy source

used to treat prostate cancer


candidate for intravascular use






Run 73


1st
insertion
(brass)




Run 74

2nd
insertion
(brass)









Run 74 calibration



Run 74 fit







Some concluding thoughts ...

- ✦ **Calorimetry is SLOW**
 - needs long time to thermally stabilize**
 - typically need multiple determinations**
 - different / absorbers / Monte Carlo calc. verifications**
- ✦ **Calorimetry NOT Necessarily Non-Destructive method**
- ✦ **NOW, accuracy is in range of ± 1 or 2 percent**
 - Largely due to baseline instabilities and**
 - uncertainties in establishing baselines to get ΔP**
- ✦ **Power may be measured very accurately**
 - But still need average energy per decay to get Activity**

END