

A Model for Analyzing Components of Uncertainty Encountered in ^3H -Standard Efficiency Tracing in $4\pi\beta$ Liquid Scintillation Counting

Brian E. Zimmerman and Ron Collé


*Physics Laboratory, National Institute of Standards and Technology
Gaithersburg, MD USA*

Winter Meeting of the American Nuclear Society
Washington, D.C., 13 November 2000

Scope of Presentation

- Goals
 - Introduce ISO philosophy of uncertainty analysis
 - Develop metrological mindset of searching for uncertainty
 - Apply ISO guidelines to development of uncertainty model for LS efficiency tracing techniques
- What is not covered
 - How to analyze all the components, experiment design, etc.
 - What to do with correlations
- Example: Standardization of ^{63}Ni (NIST SRM 4226C)

What is a *Measurement Model*?



“Formalization of the relationship between input and output quantities involved in a measurement and how the uncertainty in an input quantity leads to an uncertainty in the output quantity”


ISO Recommendations

- Determine mathematical expression for relationship between measurand Y , from individual input quantities X_i .
 - The input quantities may include variables or uncertainty components that are not explicitly needed to calculate the estimated measurand, y .
 - Example: $R = A \varepsilon d \cdot (xyz...)$
- Determine estimated value of estimated input quantities, x_i , and “standard uncertainties”, u_i .
 - u_i are expressed as estimated standard deviations (or standard deviations of the mean).
 - Uncertainty components assumed to correspond to standard deviations, irrespective of evaluation method. (Type A and B).

ISO Recommendations, Cont'd

- Evaluate covariances.
- Calculate y and “combined standard uncertainty”, u_c (u_i combined in quadrature)
- Form “expanded uncertainty”, $U = k \cdot u_c$, where k is coverage factor. Limits are supposed to correspond to a particular confidence level.
- Report y , u_c (*or* U), specification of value and assumptions for choice of k , and complete specification of all uncertainty components, u_i .

Simple Case



Component Analysis

- Just because you can think of an effect doesn't mean that the effect is present
 - Statistical analysis required
 - ANOVA can help
- An effect may be present in one experiment and not in another
- Effects may be embodied in other components

Example: LS Composition Effects

- Prepare two sets of six cocktails, each set having different cocktail composition (aqueous fraction, carrier concentration, etc.)
- Repeat for second radionuclide
- Compare within-measurement and between-measurement variabilities to determine presence of effect.
- ANOVA

Example: Test for Composition Effect

Nuclide X

Trial	$C_A, \text{Comp. A}$	$C_A, \text{Comp. B}$
1	122.6	124.1
2	123.4	124.5
3	123.1	124.6
4	122.8	124.3
5	123.3	124.7
6	123.2	124.2
Mean	123.1	124.4
s	0.3 (0.25 %)	0.2 (0.19 %)

Nuclide Y

Trial	$C_A, \text{Comp. A}$	$C_A, \text{Comp. B}$
1	523.1	523.4
2	520.3	519.8
3	521.6	522.1
4	522.1	523.3
5	522.4	521.2
6	520.5	523.5
Mean	521.7	522.2
s	1.1 (0.21 %)	1.5 (0.28 %)

ANOVA Analysis

Nuclide X

122.6	124.1		Anova: Single Factor					
123.4	124.5							
123.1	124.6		SUMMARY					
122.8	124.3		<i>Groups</i>	<i>Count</i>	<i>Sum</i>	<i>Average</i>	<i>Variance</i>	
123.3	124.7		Column 1	6	738.4	123.0667	0.094667	
123.2	124.2		Column 2	6	746.4	124.4	0.056	
			ANOVA					
			<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>
			Between Groups	5.333333	1	5.333333	70.79646	7.55E-06
			Within Groups	0.753333	10	0.075333		4.964591
			Total	6.086667	11			

Composition effect present (95 % confidence)


ANOVA Analysis, cont'd

Nuclide Y

523.1	523.4		Anova: Single Factor						
520.3	519.8								
521.6	522.1		SUMMARY						
522.1	523.3		<i>Groups</i>	<i>Count</i>	<i>Sum</i>	<i>Average</i>	<i>Variance</i>		
522.4	521.2		Column 1	6	3130	521.6667	1.202667		
520.5	523.5		Column 2	6	3133.3	522.2167	2.221667		
			ANOVA						
			<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
			Between Groups	0.9075	1	0.9075	0.53003	0.483293	4.964591
			Within Groups	17.12167	10	1.712167			
			Total	18.02917	11				

No composition effect present (95 % confidence)

Two-Phototube Coincidence Liquid Scintillation Counting with the CIEMAT/NIST Efficiency Tracing Method




$$\varepsilon = \int_0^{E_{\max}} \left\{ 1 - \exp \left[-M^{-1}EQ(E)W(E) \right] \right\}^2 P(Z,E) dE \times \left(\int_0^{E_{\max}} P(Z,E) dE \right)^{-1}$$

$$\varepsilon = f(M)$$


CIEMAT/NIST Efficiency Tracing Method

- Prepare chemically-matched LS cocktails of ^3H and nuclide of interest over some quenching range.
- Experimental ^3H efficiency is used to determine a “Figure of Merit”, M , which is used to calculate the efficiency of nuclide of interest.


CIEMAT/NIST Efficiency Tracing Method



^3H efficiency and quench determination



Quench and rate determination for nuclide of interest



EFFY calculation

- S_i -- for n_m measurements on each ^3H , ^{63}Ni and Blank sample
 $i = 1, 2, 3, \dots, n_m$
 S_j -- for n_s samples in each efficiency tracing series
 $j = 1, 2, 3, \dots, n_s$
 S_k -- for n_e tracing series of samples
 $k = 1, 2, 3, \dots, n_e$

Conclusion

- We have developed a model for analyzing uncertainties encountered using the CIEMAT/NIST efficiency tracing method.
- Model is consistent with ISO Guide to Expression of Uncertainty in Measurement.
- Present model neglects correlated effects
- Model can evolve as more data are collected.